

Excerpt from *Mol Cell Proteomics*. 2013 Feb 6

Site specific glycoforms of haptoglobin in liver cirrhosis and hepatocellular carcinoma

Pompach P, Brnakova Z, Sanda M, Wu J, Edwards N, and Goldman R. Georgetown University, United States.

SUMMARY

Haptoglobin is a liver secreted glycoprotein with four N-glycosylation sites. Its glycosylation was reported to change in several cancer diseases which prompted us to examine site specific glycoforms of haptoglobin in liver cirrhosis and hepatocellular carcinoma. To this end, we have used two dimensional separation composed of hydrophilic interaction and nano reverse phase chromatography coupled to QTOF mass spectrometry of the enriched glycopeptides. Our results show increased fucosylation of haptoglobin in liver disease with up to six fucoses associated with specific glycoforms of one glycopeptide. Structural analysis using exoglycosidase treatment and MALDI-MS/MS of detached permethylated glycans led to the identification of Lewis Y type structures observed particularly in the pooled hepatocellular carcinoma sample. To confirm the increase of the Lewis Y structures observed by LC-MS, we have used immunoaffinity detection with Lewis Y specific antibodies. The presence of multiply fucosylated Lewis Y glycoforms of haptoglobin in the disease context could have important functional implications.

Protein digestion and exoglycosidase treatment

Hp was digested as described previously (29). Briefly, 2.5 µg of isolated Hp was resuspended in 20 µL of 50 mM NH₄HCO₃ at pH 7.8 (Sigma-Aldrich, St. Louis, MO) with 0.05% RapiGest (Waters, Milford, MA), reduced with 5 mM DTT and alkylated with 15 mM iodoacetamide (Sigma-Aldrich, St. Louis, MO). Tryptic digest (2.5 ng/µL) (Promega, Madison, WI) was carried out at 37°C in Barocycler NEP2320 (Pressure BioSciences, South Easton, MA) or with 8 endoproteinase GluC (60 ng/µL) (Roche, Indianapolis, IN) at 25 °C overnight. The digests were desalted on a Micro Trap peptide cartridge (Michrom Bioresources, Auburn, CA), and washed 3 times with 250 µL of 0.1% aqueous TFA (Sigma-Aldrich, St. Louis, MO). The peptides were eluted with 100 µL of 60% AcN with 0.1%TFA and the eluate was dried using a SpeedVac concentrator.

[Click on the Title for the Complete Publication](#)

[Continued on Page 3, Column 1](#)

Excerpt from *Mol Cell Proteomics*. 2013 Feb 6

Quantitative LC-MS-MRM analysis of site-specific glycoforms of haptoglobin in liver disease

Sanda M, Pompach P, Brnakova Z, Wu J, Makambi K, and Goldman R. Georgetown University, United States

SUMMARY

Development of liver disease is associated with the appearance of multiply fucosylated glycoforms of haptoglobin. To analyze the disease-related haptoglobin glycoforms in liver cirrhosis and hepatocellular carcinoma, we have optimized an LC-MS-MRM workflow for glycopeptide quantification. The final quantitative analysis included twenty four site-specific glycoforms generated by treatment of a tryptic digest of haptoglobin with $\alpha(2-3,6,8)$ neuraminidase and $\beta(1-4)$ galactosidase. The combination of LC-MS-MRM with exoglycosidase digests allowed resolution of isobaric glycoforms of the Hp-T3 glycopeptide for quantification of the multiply fucosylated Lewis Y-containing glycoforms we have identified in the context of liver disease. Fourteen multiply fucosylated of the twenty examined fucosylated glycoforms increase significantly in the liver disease group compared to healthy controls with an average 5-fold increase in intensity ($p < 0.05$). At the same time, two tri-antennary glycoforms without fucoses do not increase in the liver disease group and two tetra-antennary glycoforms without fucoses show a marginal at most 40% increase in intensity. Our analysis of thirty individual patient samples (10 healthy controls (CTL), 10 cirrhosis (CIR) and 10 hepatocellular carcinoma (HCC) patients) shows that these glycoforms are substantially increased in a small subgroup of liver disease patients but do not significantly differ between the groups of HCC and CIR patients. The tri- and tetra-antennary singly fucosylated glycoforms are associated with MELD score and low platelet counts ($p < 0.05$). The exoglycosidase-assisted LC-MS-MRM workflow, optimized for the quantification of fucosylated glycoforms of haptoglobin, can be used for quantification of these glycoforms on other glycopeptides with appropriate analytical behavior.

[Click on the Title for the Complete Publication](#)

[Continued on Page 3, Column 2](#)

SURVEY 1: ACCELERATED PROTEIN DIGESTION

We are conducting a short survey to research the needs of scientists who would benefit from an innovative proteolytic enzyme (e.g. trypsin) digestion process that dramatically reduces the time needed for digestion and that works particularly well for samples that are difficult to digest and/or sequence.

The survey should take less than 5 minutes to complete and participants who complete the survey will be entered for a chance to win a \$100 Amazon gift certificate.

Thanks for your participation and good luck!

<https://www.surveymonkey.com/s/trypsin-digestion-survey>

[Site specific glycoforms of haptoglobin in liver cirrhosis and hepatocellular carcinoma](#): Continued from Page 2, Column 1

Determination of Glyco-site occupancy

Occupancies of all four Hp glycosites were quantified by comparison of XIC precursor ion intensities of unoccupied N-(asparagine) peptides and deglycosylated D-(aspartic acid) peptides following PNGaseF deglycosylation under 180 water as described previously (30, 31). Briefly, Hp (2 pmols) was digested after reduction and alkylation by trypsin or endoproteinase GluC as described above. The GluC digest was split in half and one portion was further digested with trypsin while the other was processed as described below; this is done to generate singly glycosylated proteolytic peptides of all four Hp glycosites. All enzymatic digests (tryptic, GluC and combination) were heated for 15 mins at 90°C to deactivate proteolytic enzymes and 11 evaporated to dryness using speedvac. Separately, 2µl of G7 reaction buffer and 0.5 µl of PNGase F (New England BioLabs, Ipswich, MA) were evaporated and diluted in 20 µl of 180water (Cambridge isotope Laboratories, Andover, MA). This solution was used to dissolve the dried peptides for de-glycosylation (1.5 hour at 37°C) using Barocycler NEP2320 (Pressure BioSciences, South Easton, MA). De-glycosylated 180 labeled peptides (2pmols) were analyzed directly after the PNGaseF enzymatic treatment on a QTOF mass analyzer using an IDA workflow.

DISCUSSION

In summary, we have detected multiply fucosylated glycoforms of Hp with up to six fucoses associated with the presence of Le(y)-type glycoforms of the T3 glycopeptide in the HCC disease group. Analysis of detached permethylated N-glycans and immunostaining with anti-Le(y) antibody confirmed the presence of this glycoform. Function of these unusual minor glycoforms is at present unknown; further studies would be required to establish association of the glycoforms with viral infection or liver disease. The hyper fucosylation likely affects interactions with lectins and other binding partners which could impact the disease progression of viral infection to HCC.

[Quantitative LC-MS-MRM analysis of site-specific glycoforms of haptoglobin in liver disease](#): Continued from Page 2, Column 2

Quantification of Hp glycopeptides

Hp isolated from patient samples (2.5 µg based on peak area of the RP HPLC chromatogram) was spiked with labeled internal standard (VVLHPNYSQVDIGLIK), reduced with DTT, alkylated with IAA, dissolved in 50mM ammonium bicarbonate buffer (Sigma-Aldrich, St. Louis, MO) containing 0.05% RapiGest, and digested with trypsin (Promega, Madison, WI) using Barocycler NEP2320 (Pressure BioSciences, South Easton, MA) at 37°C. Tryptic digests were first treated with α(2-3,6,8)neuraminidase (New England BioLabs, Ipswich, MA) at 37°C in 50mM sodium acetate, 5mm CaCl₂, pH 5.5 (New England BioLabs, Ipswich, MA). When appropriate, the digest with α(2-3,6,8)neuraminidase was followed by β(1-4)galactosidase (New England BioLabs, Ipswich, MA) overnight at 37°C in 50mM sodium acetate, 5mm CaCl₂, pH 5.5 (New England BioLabs, Ipswich, MA).

CONCLUSION

We have developed an exoglycosidase-assisted LC-MS-MRM workflow for quantification of glycoforms of Hp-T3 in liver disease context. The incorporation of exoglycosidase digests into the LC-MS-MRM workflow allows separation of isobaric fucosylated glycoforms associated with cirrhosis and HCC. The fucosylated glycoforms of Hp-T3 are up-regulated on average 5-28 fold in liver disease but the Lewis Y-type glycoforms increase in perhaps 10-20% of the patients. This is not sufficient for a consistent detection of the patients with HCC on the background of liver cirrhosis. Additional appropriately designed studies are needed to determine whether this subset of patients has different prognosis or shares common clinical characteristics. Glycopeptide MRM assays require analyte-specific optimization. To facilitate this task, we suggest a workflow, including fragment-specific CE calculations, for optimization of the LC-MS-MRM analysis of complex N-glycopeptides. Applicability of the workflow is at present limited to isolated N-glycoproteins or their simple mixtures to allow resolution of glycoforms monitored by the common oxonium ions. While general applicability of the assays has limitations, we expect that the exoglycosidase-assisted LC-MS-MRM will provide an excellent tool for verification studies comparing site specific glycoforms of interest in relevant disease conditions.

The Barocycler NEP2320



[Click Here for PCT Applications](#)

SURVEY 2: CELL DISRUPTION

We are conducting survey to research the needs of scientists who would benefit from an innovative cell disruption technology platform that provides manual and fully-automated cell disruption from milliliter to liter volumes.

The survey should take less than 5 minutes to complete and participants who complete the survey will be entered for a chance to win a \$100 Amazon gift certificate.

Thank you for participating and good luck!

<http://www.surveymonkey.com/s/cell-disruption-survey>

Presented at the Mid America Forensic DNA Conference,
April, 2013

**Pressure Cycling Technology (PCT) Reduces
Effects of Inhibitors of the PCR**

Pam Marshall, University of North Texas Health Science Center

PCT

- While investigating the use of PCT to attempt to increase DNA yield in challenged samples, we observed that PCT reduced the effects of inhibition on downstream DNA analyses.
- I hypothesized that conditions of extreme pressure alter the conformation of some inhibitors, rendering them less effective at reducing the yield of PCR products in forensic DNA analysis.



**Pressure Cycling Equipment
and Consumables**



- Barocycler@ NEP3229**
- Bench top instrument
 - Capable of processing up to three samples simultaneously using PULSE Tubes or up to 48 samples using MicroTubes



- Disc-less PULSE Tubes**
- Pressure Used to Lyse Samples for Extraction hold up to 1.2 mL
 - During the PCT process, the Ram compresses the sample and processing buffer against the inside of the Cap, resulting in pressure being transmitted to the sample



- PCT MicroTubes with Cartridge Adapter**
- Specially designed single use tubes are available in three sizes (50, 100, and 150µl)




DNA INHIBITORS

Two potent PCR inhibitors, hematin and humic acid, were evaluated

Hematin


- Metal chelating molecule found in RBCs
- Forms a stable complex with the DNA polymerase
- May also cause a dissociation of the DNA-polymerase complex, thereby inhibiting polymerase activity



www.chemicalbook.com

Humic acid

- Component in soil
- Often encountered in samples that have been buried, such as skeletal remains
- Proposed mechanisms: chelates with magnesium ions needed for DNA polymerase or inhibits the PCR via sequence specific binding to DNA, thus limiting the amount of available template



Click on the Title for the Complete Presentation

Continued on Page 4, Column 2

Pressure Cycling Technology (PCT) Reduces Effects of Inhibitors of the PCR: Continued from Page 4, Column 1

**NPC vs PCT:
Number of Alleles Detected**

	Sample Name	Total RFUs in Sample	Average PH w/n Sample	Standard Deviation PH w/n Sample	# of Alleles Detected
Biom 1	NPC 2.001.1	2318	208.8698887	129.3208368	12
	PCT 2.001.1	2784	120.1904782	51.20211669	21
	NPC 2A.001.1	2460	164	144.5422294	15
	PCT 2A.001.1	2608	141.0284615	75.80643252	26
Biom 2	NPC 2.002.1	11137	445.48	121.8948925	25
	PCT 2.002.1	18328	853.04	235.1239175	25
	NPC 2A.002.1	7005	280.2	147.6211367	25
	PCT 2A.002.1	8124	244.95	178.8058774	25
Biom 3	NPC 2.003.1	736	92	58.81690714	8
	PCT 2.003.1	1927	81.28571429	38.28987159	21
	NPC 2A.003.1	612	102	76.57853588	6



Discussion

- The results support that pressure cycling technology (PCT) reduces inhibitory effects and thus enhances yield of amplified products of both hematin and humic acid inhibited samples.
- Based on the results obtained in this study, this method can improve the recovery of alleles from challenged or inhibited DNA samples.



Conclusions

- Results demonstrate enhanced PCR efficiency for inhibited samples following PCT compared with those samples not exposed to pressure technology.
- These results indicate as a proof of concept that PCT is a viable method to overcome the inhibitory effects of hematin and HA on the PCR.
- This research study demonstrates the capabilities and potential of PCT applications for forensic DNA analysis of biological evidentiary samples.

